

| Categoria Principale | Categoria Secondaria | N° Reg. / CND / N° Repertorio | Nome commerciale | Azienda | Misure cmxcm (n°pz x conf.) | Cod. Rimb. |
|----------------------|----------------------|-------------------------------|---------------------|---------------|-----------------------------|------------|
| Wound Care | Film Poliuretano | M04040601 | Protect Film | Farmac-Zabban | 4x5 (100) | -- |
| Wound Care | Film Poliuretano | M04040601 | Protect Film | Farmac-Zabban | 6x7(100) | -- |
| Wound Care | Film Poliuretano | M04040601 | Protect Film | Farmac-Zabban | 10x12(100) | -- |
| Wound Care | Film Poliuretano | M04040601 | Protect Film | Farmac-Zabban | 10x25(50) | -- |
| Wound Care | Film Poliuretano | M04040601 | Protect Film | Farmac-Zabban | 15x20(50) | -- |
| Wound Care | Film Poliuretano | M04040601 | Protect Film | Farmac-Zabban | 20x30(20) | -- |
| Wound Care | Film Poliuretano | M04040601 | Protect Film Roll | Farmac-Zabban | 5x10m(1) | -- |
| Wound Care | Film Poliuretano | M04040601 | Protect Film Roll | Farmac-Zabban | 10x10m(1) | -- |
| Wound Care | Film Poliuretano | M04040601 | Protect Film Roll | Farmac-Zabban | 10x15m(1) | -- |
| Wound Care | Film Poliuretano | M04040601 | Euroderm | Eurofarm | 4x5(100) | |
| Wound Care | Film Poliuretano | M04040601 | Euroderm | Eurofarm | 6x7(100) | |
| Wound Care | Film Poliuretano | M04040601 | Euroderm | Eurofarm | 10x12(50) | |
| Wound Care | Film Poliuretano | M04040601 | Euroderm | Eurofarm | 15x20(10) | |
| Wound Care | Film Poliuretano | M04040601 | Euroderm | Eurofarm | 10x25(10) | |
| Wound Care | Film Poliuretano | M04040601 | Euroderm Roll | Eurofarm | 2,5x10m (408) | |
| Wound Care | Film Poliuretano | M04040601 | Euroderm Roll | Eurofarm | 5x10m(204) | |
| Wound Care | Film Poliuretano | M04040601 | Euroderm Roll | Eurofarm | 10x10m(102) | |
| Wound Care | Film Poliuretano | M04040601 | Euroderm Roll | Eurofarm | 15x10m(72) | |
| Wound Care | Film Poliuretano | M04040601 | Euroderm Roll | Eurofarm | 20x10m(52) | |
| Wound Care | Film Poliuretano | 903676322 M040406 | Bioclusive | Acelity | 5 x 7,5 (100) | -- |
| Wound Care | Film Poliuretano | 903595573 M040406 | Bioclusive | Acelity | 10 x 12,5 (50) | -- |
| Wound Care | Film Poliuretano | 903676334 M040406 | Bioclusive | Acelity | 12,5 x 17,5 (20) | -- |
| Wound Care | Film Poliuretano | 903676385 M040406 | Bioclusive | Acelity | 10 x 25 (20) | -- |
| Wound Care | Film Poliuretano | 903676447 M040406 | Bioclusive | Acelity | 20 x 25 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 | Film Polyskin II | Covidien | 5 x 7 (100) | -- |
| Wound Care | Film Poliuretano | M04040601 | Film Polyskin II | Covidien | 10 x 12 (50) | -- |
| Wound Care | Film Poliuretano | M04040601 | Film Polyskin II | Covidien | 15 x 20 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 | Film Polyskin II | Covidien | 10 x 20 (20) | -- |
| Wound Care | Film Poliuretano | M04040601 | Film Polyskin II | Covidien | 20 x 25 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 | Film Polyskin II | Covidien | 4 x 4 (100) | -- |
| Wound Care | Film Poliuretano | M04010202 | Fixomull Trasparent | BSN | 5 x 10 m (1) | -- |
| Wound Care | Film Poliuretano | M04010202 | Fixomull Trasparent | BSN | 10 x 10 m (1) | -- |
| Wound Care | Film Poliuretano | M04010202 | Fixomull Trasparent | BSN | 10 x 2 m (1) | -- |
| Wound Care | Film Poliuretano | M04010202 | Fixomull Trasparent | BSN | 15 x 10 m (1) | -- |
| Wound Care | Film Poliuretano | M04040601 174529R | Hydrofilm | Hartmann | 6 x 7 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 174544R | Hydrofilm | Hartmann | 10 x 12,5 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 174546R | Hydrofilm | Hartmann | 10 x 15 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 174547R | Hydrofilm | Hartmann | 15 x 20 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 174549R | Hydrofilm | Hartmann | 20 x 30 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 | Hydrofilm | Hartmann | 10 x 25 (25) | -- |
| Wound Care | Film Poliuretano | M04040601 174542R | Hydrofilm | Hartmann | 6 x 7 (100) | -- |
| Wound Care | Film Poliuretano | M04040601 174552R | Hydrofilm | Hartmann | 10 x 12,5 (100) | -- |
| Wound Care | Film Poliuretano | M04040601 174554R | Hydrofilm | Hartmann | 10 x 15 (50) | -- |
| Wound Care | Film Poliuretano | M04040601 174556R | Hydrofilm | Hartmann | 15 x 20 (50) | -- |

| | | | | | | |
|------------|------------------|----------------------|-------------------|-----------|----------------|----|
| Wound Care | Film Poliuretano | M04040601 174559R | Hydrofilm | Hartmann | 10 x 25 (25) | -- |
| Wound Care | Film Poliuretano | M04040601 174561R | Hydrofilm | Hartmann | 12 x 25 (25) | -- |
| Wound Care | Film Poliuretano | M04010102 175364R | Hydrofilm plus | Hartmann | 5 x 7,2 (5) | -- |
| Wound Care | Film Poliuretano | M04010102 175403R | Hydrofilm plus | Hartmann | 9 x 10 (5) | -- |
| Wound Care | Film Poliuretano | M04010102 175405R | Hydrofilm plus | Hartmann | 9 x 15 (5) | -- |
| Wound Care | Film Poliuretano | M04010102 175407R | Hydrofilm plus | Hartmann | 10 x 20 (5) | -- |
| Wound Care | Film Poliuretano | M04010102 175491R | Hydrofilm plus | Hartmann | 5 x 7,2 (50) | -- |
| Wound Care | Film Poliuretano | M04010102 175495R | Hydrofilm plus | Hartmann | 9 x 10 (50) | -- |
| Wound Care | Film Poliuretano | M04010102 175500R | Hydrofilm plus | Hartmann | 9 x 15 (25) | -- |
| Wound Care | Film Poliuretano | M04010102 175503R | Hydrofilm plus | Hartmann | 10 x 12 (25) | -- |
| Wound Care | Film Poliuretano | M04010102 175507R | Hydrofilm plus | Hartmann | 10 x 20 (25) | -- |
| Wound Care | Film Poliuretano | M04010102 175511R | Hydrofilm plus | Hartmann | 10 x 25 (25) | -- |
| Wound Care | Film Poliuretano | M04010102 175516R | Hydrofilm plus | Hartmann | 10 x 30 (25) | -- |
| Wound Care | Film Poliuretano | M04010102 178766R | Hydrofilm Roll | Hartmann | 5 x 10m (1) | -- |
| Wound Care | Film Poliuretano | M04010102 178799R | Hydrofilm Roll | Hartmann | 10 x 2m (1) | -- |
| Wound Care | Film Poliuretano | M04010102 178802R | Hydrofilm Roll | Hartmann | 10 x 10m (1) | -- |
| Wound Care | Film Poliuretano | M04010102 178803R | Hydrofilm Roll | Hartmann | 15 x 10m (1) | -- |
| Wound Care | Film Poliuretano | M04040601 | Leukomed T | BSN | 7,2 x 5 (50) | -- |
| Wound Care | Film Poliuretano | M04040601 | Leukomed T | BSN | 8 x 10 (50) | -- |
| Wound Care | Film Poliuretano | M04040601 | Leukomed T | BSN | 11 x 14 (50) | -- |
| Wound Care | Film Poliuretano | M04040601 | Leukomed T | BSN | 10 x 12 (50) | -- |
| Wound Care | Film Poliuretano | M04040601 | Leukomed T | BSN | 10 x 25 (50) | -- |
| Wound Care | Film Poliuretano | M04040601 | Leukomed T | BSN | 15 x 20 (50) | -- |
| Wound Care | Film Poliuretano | M04040601 | Leukomed T | BSN | 15 x 25 (50) | -- |
| Wound Care | Film Poliuretano | M040499 506435/R | Mepitel Film | Molnlycke | 6,5 x 7 (10) | |
| Wound Care | Film Poliuretano | M040499 1018238/R | Mepitel Film | Molnlycke | 10,5 x 12 (10) | |
| Wound Care | Film Poliuretano | M040499 1018239/R | Mepitel Film | Molnlycke | 10,5 x 25 (10) | |
| Wound Care | Film Poliuretano | M040499 1018240/R | Mepitel Film | Molnlycke | 15,5 x 20 (10) | |
| Wound Care | Film Poliuretano | M04010202 156872R | Mepore Film | Molnlycke | 6 x 7 (100) | -- |
| Wound Care | Film Poliuretano | M04010202 1022864/RR | Mepore Film | Molnlycke | 10 x 12 (70) | -- |
| Wound Care | Film Poliuretano | M04010202 1025865/R | Mepore Film | Molnlycke | 10 x 25 (20) | -- |
| Wound Care | Film Poliuretano | M04010202 1022866/R | Mepore Film | Molnlycke | 15 x 20 (10) | -- |
| Wound Care | Film Poliuretano | M04010202 1022867/R | Mepore Film | Molnlycke | 20 x 30 (5) | -- |
| Wound Care | Film Poliuretano | M04010202 | Mepore Film | Molnlycke | 6 x 7 (10) | -- |
| Wound Care | Film Poliuretano | M04010202 | Mepore Film | Molnlycke | 10 x 12(10) | -- |
| Wound Care | Film Poliuretano | M04010202 | Mepore Film | Molnlycke | 10 x 25 (10) | -- |
| Wound Care | Film Poliuretano | M04010102 156240R | Mepore Film & Pad | Molnlycke | 4 x 5 (85) | -- |
| Wound Care | Film Poliuretano | M04010102 | Mepore Film & Pad | Molnlycke | 4 x 5 (5) | -- |
| Wound Care | Film Poliuretano | M04010102 1023855/R | Mepore Film & Pad | Molnlycke | 5 x 7 (85) | -- |
| Wound Care | Film Poliuretano | | Mepore Film & Pad | Molnlycke | 5 x 7 (5) | |
| Wound Care | Film Poliuretano | M04010102 1023856/R | Mepore Film & Pad | Molnlycke | 9 x 10 (30) | -- |
| Wound Care | Film Poliuretano | M04010102 | Mepore Film & Pad | Molnlycke | 9 x 10 (5) | -- |
| Wound Care | Film Poliuretano | M04010102 1023857/R | Mepore Film & Pad | Molnlycke | 9 x 15 (30) | -- |
| Wound Care | Film Poliuretano | M04010102 | Mepore Film & Pad | Molnlycke | 9 x 15 (5) | -- |
| Wound Care | Film Poliuretano | M04010102 1023858/R | Mepore Film & Pad | Molnlycke | 9 x 20 (30) | -- |

| | | | | | | |
|------------|------------------|----------------------|-------------------------|--------------------|---------------|--------------|
| Wound Care | Film Poliuretano | M04010102 | Mepore Film & Pad | Molnlycke | 9 x 20 (5) | -- |
| Wound Care | Film Poliuretano | M04010102 1023859/R | Mepore Film & Pad | Molnlycke | 9 x 25 (30) | -- |
| Wound Care | Film Poliuretano | M04010102 1023860/R | Mepore Film & Pad | Molnlycke | 9 x 30 (25) | -- |
| Wound Care | Film Poliuretano | M04010102 102361/R | Mepore Film & Pad | Molnlycke | 9 x 35 (25) | -- |
| Wound Care | Film Poliuretano | M04010102 156240R | Mepore Film & Pad Ovale | Molnlycke | 5 x 7 (85) | -- |
| Wound Care | Film Poliuretano | M04010102 | Mepore Film & Pad Ovale | Molnlycke | 5 x 7 (5) | -- |
| Wound Care | Film Poliuretano | M040499 | Opsite Flexifix | Smith&Nephew | 5 x 1000 (1) | -- |
| Wound Care | Film Poliuretano | 906436769 M040499 | Opsite Flexifix | Smith&Nephew | 10 x 1000 (1) | -- |
| Wound Care | Film Poliuretano | 906436757 M040499 | Opsite Flexifix | Smith&Nephew | 15 x 1000 (1) | -- |
| Wound Care | Film Poliuretano | 901603858 M040499 | Opsite Flexigrid | Smith&Nephew | 6 x 7 (100) | -- |
| Wound Care | Film Poliuretano | 909118871 M040499 | Opsite Flexigrid | Smith&Nephew | 10 x 12 (10) | -- |
| Wound Care | Film Poliuretano | M040499 | Opsite Flexigrid | Smith&Nephew | 10 x 12 (50) | -- |
| Wound Care | Film Poliuretano | M040499 | Opsite Flexigrid | Smith&Nephew | 15 x 20 (10) | -- |
| Wound Care | Film Poliuretano | M040499 | Opsite Flexigrid | Smith&Nephew | 12 x 25 (20) | -- |
| Wound Care | Film Poliuretano | M04040601 57281R | Suprasorb F | Smith&Nephew | 5 x 7 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 57281R | Suprasorb F | Smith&Nephew | 5 x 7 (100) | -- |
| Wound Care | Film Poliuretano | M04040601 57281R | Suprasorb F | Smith&Nephew | 10 x 12 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 57281R | Suprasorb F | Smith&Nephew | 10 x 12 (50) | -- |
| Wound Care | Film Poliuretano | M04040601 57281R | Suprasorb F | Lohmann & Rauscher | 10 x 25 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 96788R | Suprasorb F | Lohmann & Rauscher | 15 x 20 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 96788R | Suprasorb F | Lohmann & Rauscher | 20 x 30 (10) | -- |
| Wound Care | Film Poliuretano | M04040601 57372R | Suprasorb F in rotolo | Lohmann & Rauscher | 10 x 1m (1) | -- |
| Wound Care | Film Poliuretano | M04040601 57372R | Suprasorb F in rotolo | Lohmann & Rauscher | 10 x 10m (1) | -- |
| Wound Care | Film Poliuretano | M04040601 57372R | Suprasorb F in rotolo | Lohmann & Rauscher | 15 x 10m (1) | -- |
| Wound Care | Film Poliuretano | 23093R M04040601 | Tegaderm film | 3M | 4 x 4 (100) | -- |
| Wound Care | Film Poliuretano | 23083R M04040601 | Tegaderm film | 3M | 6 x 7 (100) | -- |
| Wound Care | Film Poliuretano | 23093R M04040601 | Tegaderm film | 3M | 10 x 12 (50) | 09.21.12.003 |
| Wound Care | Film Poliuretano | 23093R M04040601 | Tegaderm film | 3M | 10 x 25 (20) | -- |
| Wound Care | Film Poliuretano | 23093R M04040601 | Tegaderm film | 3M | 15 x 20 (10) | 09.21.12.006 |
| Wound Care | Film Poliuretano | 23093R M04040601 | Tegaderm film | 3M | 20 x 30 (10) | -- |
| Wound Care | Film Poliuretano | 15738R M0599 | Tegaderm Roll | 3M | 5 x 10 m (4) | -- |
| Wound Care | Film Poliuretano | 15738R M0599 | Tegaderm Roll | 3M | 10 x 10 m (4) | -- |
| Wound Care | Film Poliuretano | 15738R M0599 | Tegaderm Roll | 3M | 15 x 10 m (4) | -- |
| Wound Care | Film Poliuretano | M04010202 347822/R | Mepore Film Roll | Molnlycke | 5 x 10m (1) | |
| Wound Care | Film Poliuretano | M04010202 1018221/R | Mepore Film Roll | Molnlycke | 10 x 10m (1) | |
| Wound Care | Film Poliuretano | M04010202 1018222/R | Mepore Film Roll | Molnlycke | 15 x 10m (1) | |
| Wound Care | Film Poliuretano | M0599 15738/R | Tegaderm Film | 3M | 4x4(100) | -- |
| Wound Care | Film Poliuretano | M04040601 23083 | Tegaderm Film | 3M | 6x7(100) | -- |
| Wound Care | Film Poliuretano | M04040601 23093 | Tegaderm Film | 3M | 10x12(50) | -- |
| Wound Care | Film Poliuretano | M04040601 23093 | Tegaderm Film | 3M | 10x25(20) | -- |
| Wound Care | Film Poliuretano | M04040601 23093 | Tegaderm Film | 3M | 15x20(10) | -- |
| Wound Care | Film Poliuretano | M04040601 23093 | Tegaderm Film | 3M | 20x30(10) | -- |
| Wound Care | Film Poliuretano | M04040601 23093 | Tegaderm Film | 3M | 5x10mt(4) | -- |
| Wound Care | Film Poliuretano | M04040601 23093 | Tegaderm Film | 3M | 10x10mt(4) | -- |
| Wound Care | Film Poliuretano | M04040601 23093 | Tegaderm Film | 3M | 15x10mt(4) | -- |